

COURIERS

Passive temperature management courier solutions with VIP insulation



Technical specifications

The CSafe AcuTemp® courier solutions, with Vacuum Insulated Panel (VIP) technology, are designed for physicians, pharmacists, researchers, pharmaceutical representatives and other healthcare professionals making “last-mile” delivery of temperature-sensitive life science products. The AcuTemp couriers are simplistic, customizable and reusable temperature management solutions that offer options to accommodate a variety of payload sizes, shipping durations and temperature ranges.

Peace of Mind® thermal protection for vital healthcare product deliveries

Product features:

- Superior ThermoCor® VIP insulation, providing up to 10 times greater thermal efficacy per inch than standard foam insulation
- Designs to meet four different temperature ranges:
 - ▶ Frozen (-20°C)
 - ▶ +2°C to +8°C
 - ▶ +2°C to +25°C
 - ▶ +15°C to +25°C
- Broad portfolio available in three sizes: 3L, 6L and 12L
- Simplistic conditioning and pack-out design, with color-coded refrigerant belts
- Lightweight, highly portable and easy to clean
- Durable hard rubber handles on top, front and sides for convenience and flexible portability
- Water and skid resistant bottom
- Side and front pockets for convenient storage
- Adjustable and removable shoulder straps for comfortable hand or shoulder carrying (available on 12L)
- Configurable with most wheeled carts using secure bottom strap (available on 12L)

Safeguarding vital payloads for last-mile delivery

- Ultimate temperature protection of pharmaceuticals during final stages of transit
- Reliable temperature control across four different temperature ranges
- Efficient and unique refrigerant design allowing for pack-out simplicity
- Time-saving solution without any ‘bench time’ required for refrigerant conditioning




AcuTemp Courier solutions

Technical specifications (continued)

Unmatched thermal protection with ThermoCor VIP:

- VIP enables the couriers to effectively and efficiently maintain payload temperature
- Superior VIP insulation provides up to 10 times greater thermal efficiency per-inch than standard foam insulation
- Best-in-class performance of VIP insulation contributes to the protection of vital payloads

ThermoCor® VIP achieves far superior thermal efficiency compared to commonly used insulation



Simplistic pack-out design:



Operational specifications:

	Temperature Ranges			
	Frozen (-20°C)	+2°C to +8°C	+2°C to +25°C	+15°C to +25°C
Duration	12 hours	12 hours 48 hours	12 hours	24 hours
Profile	Universal	Seasonal Universal	Universal	Universal



Keeping temperatures precisely where they need to be

US Toll-Free: 1.888.323.9576 | International: +1.937.245.6350
Email: Sales@CSafeGlobal.com | Web: CSafeGlobal.com

CSafe, ThermoCor, AcuTemp and the logo designs presented in this material are trademarks or registered trademarks of CSafe Global. © 2019 CSafe.

DSCOURIER071619US

