



CGTULTRA

Optimal Thermal Protection

CSafe's CGT ULTRA portfolio mitigates the challenges of real-world distribution with flexible solutions designed to solve the wide range of cold chain requirements in cell and gene therapies.

Our CGT ULTRA cell and gene solutions provide industry-leading temperature protection with vacuum insulated panel (VIP) insulation and extended hold times. In addition, these solutions have been designed to hibernate in the case of unexpected delays. Simplified product handling and integrated features help avoid errors during the pack-out process making chain of custody checks easier.

With your commercialization process in mind, CSafe has designed this product portfolio to meet the rigorous Cell & Gene supply chain needs.



Key Features



Superior VIP insulation, providing optimal thermal protection



Dedicated solutions for cold and frozen payload needs



Range of payload sizes available: 5L, 8L, 16L and 35L



Extended hold times offering thermal protection for 120+hr shipment durations



Shipping delay protection



Integration of telemetry and online customer portal access



Lightweight design with integrated handles for hand-carry delivery applications



Robust, reusable solutions designed for sustainability

