



SoftboxVIP

Reusable Parcel Shippers

CSafe's portfolio of reusable parcels with vacuum insulated panel (VIP) technology, assures greater peace of mind in your cold chain.

Promising industry-leading temperature management performance and versatility for the safe transport of vital life science products, Softbox VIP shippers offer protection, simplicity and efficiency by design in pre-qualified configurations.



Key Features



High-performance reusable solution




Qualified against ISTA 7D summer and winter profiles



Best-in-class VIP insulation and PCM material for optimal thermal protection




50% longer shipment durations



Simplified pack-out process with minimal components



Eliminates need for additional packing materials



Industry-leading reuse and requalification program



Globally available via extensive supply network



PAYLOAD (LITERS)
5-56



TEMPERATURE RANGE
+2°C to +8°C
+15°C to +25°C
-25°C to -15°C
< -70°C < -80°C



DURATION (HOURS)
96-120+

PRODUCT NAME	PAYLOAD	TEMPERATURE RANGE	DURATION	INTERNAL SPACE L x W x H	EXTERNAL DIMENSIONS L x W x H	EXTERNAL DIMENSIONS (WITH TRACKER) L x W x H	TARE WEIGHT
Softbox VIP 5L	5 L 0.17 cu ft	+2°C to +8°C	96 hrs	177.8 x 177.8 x 177.8 mm 7 x 7 x 7 in	355.6 x 311.1 x 330.2 mm 14 x 12.25 x 13 in	355.6 x 311.1 x 362 mm 14 x 12.25 x 14.25 in	9.53 kg 21 lbs
			120 hrs				10.9 kg 24 lbs
		+15°C to +25°C	96 hrs				9.53 kg 21 lbs
			120 hrs				10.9 kg 22.50 lbs
		< -80°C	120 hrs				4.1 kg 9 lbs
Softbox VIP 8L	8 L 0.28 cu ft	+2°C to +8°C	96 hrs	203.2 x 203.2 x 203.2 mm 8 x 8 x 8 in	381 x 336.55 x 355.6 mm 15 x 13.25 x 14 in	381 x 336.55 x 387.4 mm 15 x 13.25 x 15.25 in	11.3 kg 25 lbs
			120 hrs				13.2 kg 29.10 lbs
		+15°C to +25°C	96 hrs				11.1 kg 24.50 lbs
			120 hrs				12.5 kg 27.50 lbs
		< -80°C	120 hrs				5 kg 11 lbs
Softbox VIP 14L	14 L 0.49 cu ft	< -70°C	240 hrs	250 x 250 x 246 mm 9.84 x 9.84 x 9.69 in	533 x 467 x 489 mm 20.98 x 18.38 x 19.52 in	n/a	14.9 kg 33 lbs
Softbox VIP 16L	16 L 0.57 cu ft	+2°C to +8°C	96 hrs	254 x 254 x 254 mm 10 x 10 x 10 in	431.8 x 387.35 x 406.4 mm 17 x 15.25 x 16 in	431.8 x 387.35 x 438.2 mm 17 x 15.25 x 17.25 in	13.8 kg 30.50 lbs
			120 hrs				15.2 kg 33.50 lbs
		+15°C to +25°C	96 hrs				13.2 kg 29.10 lbs
			120 hrs				15.2 kg 33.50 lbs
		< -80°C	120 hrs				20 kg 44 lbs
Softbox VIP 35L	35 L 1.23 cu ft	+2°C to +8°C	96 hrs	330.2 x 330.2 x 330.2 mm 13 x 13 x 13 in	508 x 463.5 x 482.6 mm 20 x 18.25 x 19 in	508 x 463.5 x 514.4 mm 20 x 18.25 x 20.25 in	20.2 kg 44.50 lbs
			120 hrs				21.5 kg 47.50 lbs
		+15°C to +25°C	96 hrs				20.2 kg 44.50 lbs
			120 hrs				21.5 kg 47.50 lbs
		< -80°C	120 hrs				27 kg 59.50 lbs
Softbox VIP 56L	56 L 1.98 cu ft	+2°C to +8°C	120 hrs	381 x 381 x 381 mm 15 x 15 x 15 in	559 x 502 x 533.5 mm 22 x 19.75 x 21 in	559 x 502 x 565.2 mm 22 x 19.75 x 22.5 in	30.4 kg 67 lbs
		-25°C to -15°C					33.1 kg 73 lbs
		< -80°C					17.3 kg 38 lbs

Contact us at passivesolutions@csafeglobal.com

CSafe has a worldwide network of comprehensive service centers, manufacturing facilities, sales offices and local hub sites to provide solutions where and when you need them. With representatives across the Americas, EMEA and APAC, local and regional support is always available.