



## Optimal Thermal Protection

The CGT Ultra provides industry-leading temperature protection with VIP insulation and extended hold times. In addition, these solutions have been designed to hibernate in the case of unexpected delays. With simplified product handling and integrated features to help avoid error during the pack-out process, your chain of custody checks are easier.



## Key Features



Superior VIP insulation, providing optimal thermal protection



Dedicated solutions for frozen payload needs



Range of payload sizes available: 2L to 56L



Extended hold times offering thermal protection for 96+hr shipment durations



Shipping delay protection and hibernation capability



Integrated tracking devices with access to CSafe Connect for visibility



Lightweight design with built-in handles for hand-carry delivery applications



Robust, reusable solutions designed for sustainability



**PAYLOAD (LITERS)**  
**2-56**



**TEMPERATURE RANGE**  
**< -60°C**  
**< -70°C**



**DURATION (HOURS)**  
**96-120**

## Technical Specifications

PRODUCT NAME	PAYLOAD	TEMPERATURE RANGE	DURATION	INTERNAL SPACE L x W x H	EXTERNAL DIMENSIONS L x W x H
CGT ULTRA 2L	2 L 0.07 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	127 x 127 x 127 mm 5 x 5 x 5 in	356 x 311 x 330 mm 14 x 12.25 x 13 in
CGT ULTRA 3L	3 L 0.13 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	152 x 152 x 152 mm 6 x 6 x 6 in	381 x 336 x 356 mm 15 x 13.25 x 14 in
CGT ULTRA 35L	35 L 1.27 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	330 x 330 x 330 mm 13 x 13 x 13 in	508 x 464 x 483 mm 22.0 x 19.75 x 22.25 in
CGT ULTRA 56L	56 L 1.95 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	381 x 381 x 381 mm 15 x 15 x 15 in	559 x 502 x 554 mm 20 x 18.25 x 19 in
CGT ULTRA 5L	5 L 0.17 cu ft	Frozen (Dry Ice) < -60°C*	120 hrs	178 x 178 x 152 mm 7 x 7 x 6 in	356 x 311 x 330 mm 14 x 12.25 x 13 in
CGT ULTRA 8L	8 L 0.30 cu ft	Frozen (Dry Ice) < -60°C*	120 hrs	203 x 203 x 203 mm 8 x 8 x 8 in	381 x 336 x 356 mm 15 x 13.25 x 14 in
CGT ULTRA 16L	16 L 0.58 cu ft	Frozen (Dry Ice) < -60°C*	120 hrs	254 x 254 x 254 mm 10 X 10 X 10 in	432 x 387 x 406 mm 17 x 15.25 x 16 in
CGT ULTRA 35L	35 L 1.27 cu ft	Frozen (Dry Ice) < -60°C*	120 hrs	330 x 330 x 330 mm 13 x 13 x 13 in	508 x 464 x 483 mm 20 x 18.25 x 19 in
CGT ULTRA 56L	56 L 1.95 cu ft	Frozen (Dry Ice) < -60°C*	120 hrs	381 x 381 x 381 mm 15 x 15 x 15 in	559 x 502 x 554 mm 22.00 x 19.75 x 22.25 in

\*-70C data available upon request

Contact us at [cellandgene@csafeglobal.com](mailto:cellandgene@csafeglobal.com)

CSafe has a global network of service centers, stations, supply points and team members to provide support and solutions where and when you need them. CSafe and the logo designs presented in this material are trademarks or registered trademarks of CSafe. ©2023 CSafe

