



Optimal Thermal Protection

The CGT Ultra provides industry-leading temperature protection with VIP insulation and extended hold times. In addition, these solutions have been designed to hibernate in the case of unexpected delays. With simplified product handling and integrated features to help avoid error during the pack-out process, your chain of custody checks are easier.



Key Features



Superior VIP insulation, providing optimal thermal protection



Dedicated solutions for frozen payload needs



Range of payload sizes available: 2L to 35L



Offering thermal protection for 96 hr shipment durations



Shipping delay protection and hibernation capability



Integrated tracking devices with access to CSafe Connect for visibility



Lightweight design with built-in handles for hand-carry delivery applications



Robust, reusable solutions designed for sustainability



PAYLOAD (LITERS)
2-35



TEMPERATURE RANGE
< -70°C



DURATION (HOURS)
96

Technical Specifications

PRODUCT NAME	PAYLOAD	TEMPERATURE RANGE	DURATION	INTERNAL SPACE L x W x H	EXTERNAL DIMENSIONS L x W x H
CGT ULTRA 2L	2 L 0.07 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	127 x 127 x 127 mm 5 x 5 x 5 in	356 x 311 x 330 mm 14 x 12.25 x 13 in
CGT ULTRA 3L	3 L 0.13 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	152 x 152 x 152 mm 6 x 6 x 6 in	381 x 336 x 356 mm 15 x 13.25 x 14 in
CGT ULTRA 35L	35 L 1.27 cu ft	Frozen (Dry Ice) < -70°C	96 hrs	330 x 330 x 330 mm 13 x 13 x 13 in	508 x 464 x 483 mm 22.0 x 19.75 x 22.25 in

*-70C data available upon request

Contact us at cellandgene@csafeglobal.com

Csafe has a global network of service centers, stations, supply points and team members to provide support and solutions where and when you need them. CSafe and the logo designs presented in this material are trademarks or registered trademarks of CSafe. ©2024 CSafe